



# UpDate!

1896 **SCHWERT SWORD** SYNONYM FÜR QUALITÄT & VERTRAUEN  
SYNONYM FOR QUALITY & TRUST



## ΟΔΟΝΤΙΑΤΡΙΚΑ-ΟΡΘΟΔΟΝΤΙΚΑ

Μικαλακοπούλου 157, Γουδή 115 27  
Τηλ: 210 7716.416, 210 7751.000, Fax: 210 7711.100  
e-mail: info@tsaprazis.gr, www.tsaprazis.gr

**SCHWERT Special Offers**

**2024**

Modern Therapy · Modern Instruments

## GRACEY curettes – all working tips are exchangeable



**3301-100 €329,00**

GRACEY Σετ, αποτελείται από:  
GRACEY 5/6, 7/8, 11/12,  
13/14, ξέστρο 3305-04S και  
δίσκο τοποθέτησης



**3301-100M5 €329,00**

GRACEY mini σετ για τις άπω επιφάνειες  
των οπίσθιων δοντιών, αποτελείται από:  
GRACEY 1/2M5, 7/8M5, 11/12M5, 13/14M5,  
ξέστρο 3305-04S και δίσκο τοποθέτησης



**3301-100M5Ti €432,00**

GRACEY σετ M5 με άκρα τιτανίου  
αποτελείται από: GRACEY 1/2M5Ti,  
7/8M5Ti, 11/12M5Ti, 13/14M5Ti  
και δίσκο τοποθέτησης



**3300-00 €44,20**

Key to exchange  
working tips


**SCHWERT**  
SWORD · ESPADA  
EPEE · SPADA

## Implant TI Curettes – exchangeable and delicate titanium tips, resharpenable and long lasting




**3302-04Ti €83,20**

COLUMBIA 4R/4L




**3303-01Ti €83,20**

LANGER L 1/2



**3303-03Ti €83,20**

LANGER L 3/4



**3303-05Ti €83,20**

LANGER L 5/6

## Composite-Instruments – ULTRA Composite



**4085-01Ti €57,20**

Εργαλείο τοποθέτησης σύνθετης  
ρητίνης, 1 FIG. 1, μικρό, 1.8 / 1.5 mm



**4085-03Ti €57,20**

Εργαλείο στίλβωσης  
σύνθετης ρητίνης, συνδυασμό  
σχήμα αχλαδιού



**4085-02Ti €57,20**

Εργαλείο τοποθέτησης  
σύνθετης ρητίνης, 1 FIG. 2, μεγάλο,  
1.8 / 1.5 mm



**4085-04Ti €57,20**

Εργαλείο συμπύκνωσης ρητίνων  
για τη διαμόρφωση μασητικών  
επιφανειών σε γομφίους, αντίθετης γωνίας.



**4085-00Ti €333,00**

Σετ εργαλείων διαμόρφωσης  
ρητίνων με τα 4 εικονιζόμενα  
εργαλεία

## LUCAS curettes



**1372-0 €74,10**

LUCAS curette, Ø 1.5 mm, non-serrated



**1372-1 €74,10**

LUCAS curette, Ø 2.5 mm, non-serrated



**1372-2 €74,10**

LUCAS curette, Ø 3.0 mm, non-serrated



**1372-10 €75,80**

LUCAS curette, Ø 1.5 mm, serrated



**1372-11 €75,80**

LUCAS curette, Ø 2.5 mm, serrated



**1372-12 €75,80**

LUCAS curette, Ø 3.0 mm, serrated

## Periosteals



**1411-05 €74,10**

BUSER, 3.8 / 4.2 mm, 17.5 cm

MADE IN GERMANY



**1411-10 €74,10**

MOLT, Fig. 3, 7.0 / 4.6 mm, 17.5 cm



**1411-24 €74,10**

MOLT, CM 4 / CM 7, 3.6 / 7.0 mm, 17.5 cm

[schwert.com](http://schwert.com)



# P-LUX instruments



617-07 €85,15

P-LUX FIG. 7, Desmotome, straight, for initial loosening of the desmodontal fibres

P-LUX Set, consisting of: 617-01, 617-03, 617-04, 617-06, 617-07, 617-20, 1372-11 and Protector-Tray 2597-08

617-100S €740,00



617-01 €85,15

P-LUX FIG. 1, straight, green, for front teeth and upper jaw



617-03 €85,15

P-LUX FIG. 3, straight, green, for front teeth and upper jaw



617-04 €85,15

P-LUX FIG. 4, curved, blue, for lower jaw



617-06 €85,15

P-LUX FIG. 6, curved, blue, for lower jaw



1. Initial loosening of the desmodontal fibres. Apply the P-LUX Desmotome first approximal then buccal and palatal. Working position pencil handle.



2. Then progressive luxation with the P-LUX instruments (first narrow then wide). Ideal positions: Approximal (mesial or distal) to the tooth which should be extracted!




3. Complete luxation with support of a wide P-LUX instrument. Gentle extraction without root fractures!



4. Inspection of the empty tooth socket. Important: Complete removal of all inflammable and connective tissue structures.

# P-LUX periotomes

with 2-axis flexibility and flexible working tip



617-20 €85,15  
straight version

617-21 €88,40  
bended version

## Cutting of Sharpey's fibres with the new SCHWERT P-LUX Periostomes.

To accomplish an atraumatic extraction the Sharpey's fibres should be separated by using a Periostome primary to the use of extraction forceps or elevators. Therefore the Periostome is driven thru the sulcus. The proper use prevents the surrounding gingiva and the periost. The SCHWERT P-LUX handle design secures a safe and tactile handling. The blade of this periostome has to be highlighted specially – **the blade is flexible and adapts perfectly to root structures, simultaneously due to its high stability it is possible to gain a certain pressure.** This is possible by choosing a special medical stainless steel and hardening treatment.

# Modified elevator BECK



615-1 €78,30  
2.8 mm

615-2 €78,30  
3.5 mm

615-3 €78,30  
4 mm





SCHWERT  
SWORD · ESPADA  
EPEE · SPADA

**i** **Deep Gripping TAC-Traction extraction forceps**

Many teeth to be extracted have already a weakened crown due to previous caries / filling treatment. Our TAC-Traction deep gripping forceps avoid complications like crown or root fractures.

**Advantages of deep gripping forceps:**

- Longer and sharpened conical jaws
- Less punctual stress on the damaged crown due to dozens of connecting points on crown & root structure
- Perfect anatomical shape for a parallel contact
- Less soft tissue trauma by subgingival placement of the jaws. Soft tissue is more retracted than squeezed.

Standard forceps

Deep gripping forceps

**140-34 €155,50**  
FIG. 34, 16.0 cm

**140-35 €155,50**  
FIG. 35, 16.0 cm

**140-67N €169,70**  
FIG. 67N, 17.5 cm

**140-36 €169,70**  
FIG. 36, 16.0 cm

**140-36A €169,70**  
FIG. 36A, 16.0 cm



**140-110**  
TAC-Traction Set consisting of deep gripping extraction forceps: **140-34, 140-35, 140-36, 140-36A, 140-67N, 140-79N** and Protector-Basket 2602-03

**€1.273,00**

**140-79N €169,70**  
FIG. 79N, 16.0 cm

### Universal root splinter forceps WITZEL

**166-300 €248,00**  
WITZEL, diamond coated, 15.5 cm

**431-10 €192,40**  
MINI-FRIEDMANN, 14.5 cm, 30°, 1.5 mm

**432 €205,00**  
LUER-FRIEDMANN, 15.5 cm, 30°, 3 mm



# SCHWERT Special Offers

2024

# UpDate!

## SuperCut scissors ULTRA SCHEREN SOISSORS


- Extremely sharp
- Less trauma

1133-400B €152,90

LA GRANGE, ULTRA, 11.5 cm

1133-400 €105,30

LA GRANGE, 11.5 cm




Standard / SuperCut



1111-400 €71,00

IRIS, straight, 11.5 cm




1117-400 €91,70

GOLDMAN-FOX, 13 cm



1125-400 €110,00

JOSEPH, 14 cm



1143-400 €100,00

KELLY, 16 cm

1111-400B €138,00

ULTRA

1117-400B €184,50

ULTRA

1125-400B €179,40

ULTRA

## Microsurgical instruments



1033-50 €57,20  
Blade holder, 15 cm



1070-00 €188,50  
Micro blade holder, 15 cm

### SAFE-EX 2 patented

– one click for more safety



1033-61 €152,00  
SAFE-EX 2 safety blade holder, 15 cm, dismantable

SCHWERT twist series

Our Needle Holder with the twisted lock avoids reliably sticking of the suture while suturing.



Suture size



2000-60 €440,00  
twist TC micro needle holder with lock, straight, 17.5 cm



2000-70 €440,00  
twist TC Kirsch, micro needle holder with lock, straight, 17.5 cm



2000-66 €451,00  
twist TC micro needle holder with lock, curved, 17.5 cm



2000-10 €439,00  
twist TC micro needle holder with lock, straight, 17.5 cm



2000-16 €439,00  
twist TC micro needle holder with lock, curved, 17.5 cm



1411-01 €74,10  
MOLT periosteal elevator, 2.6 mm / 4.0 mm, 17.5 cm



1208-300 €179,00  
Micro forceps, 1x2 teeth, straight, 17.5 cm



1208-340 €201,50  
COOLEY micro forceps, straight, 17.5 cm  
Due to atraumatic tips, perfectly suited for grasping and holding soft tissue



1135-17 €317,50  
Micro scissors, 17.5 cm, Ø 8 mm, micro serrated, curved



1137-17 €317,50  
Micro scissors, 17.5 cm, Ø 8 mm, curved, very fine

